


GeneProof VRE PCR Kit



DIRECT DETECTION FROM CLINICAL MATERIAL

- No need of previous cultivation
- Quick result from primary clinical material

EASY-TO-USE CONCEPT

- Single tube Ready-to-Use Master Mix contains all components for PCR amplification
- No additional PCR reagents pipetting necessary

SIMPLE LABORATORY WORKFLOW

- Easily combinable with other GeneProof PCR kits in one workflow

IDENTIFICATION OF *vanA* AND *vanB*

- Detects and distinguishes *vanA* and *vanB* genes
- Enables choosing appropriate therapy

CONTAMINATION PREVENTION

- Master Mix contains Uracil-DNA glycosylase (UNG) and dUTPs eliminating carryover contamination

ORDER INFORMATION

REF	PACKAGE
VRE/ISEX/025	25 reactions
VRE/ISEX/100	100 reactions



COMPATIBLE WITH A WIDE RANGE OF REAL-TIME PCR DEVICES



CERTIFIED DIAGNOSTIC TEST